



**Your Healthcare
Innovative Solution**



Your Healthcare Innovative Solution

Carelife offers comprehensive solutions for disposable medical devices, which includes product design, manufacturing, regulatory compliance.

We specialize in all phases of Private Labeling and OEM/ODM (Original Design Manufacturing). Our strict quality control follows ISO and CE standards and all of our facilities are FDA registered.

When you partner with Carelife, you receive not only a potential solution to your needs but also a global manufacturing and service team.



Choose
Carelife
Get
assurance of
Your
Life



Pen Needles

56114660162

Needle Length	Injector profile	Comments
32Gx4mm	Children Thin adults	For children and adults of BMI 20-49
31Gx5mm	Children and adolescents Thin adults Normal weight adults	Adults who use a lifted skin fold should use the 5mm under supervisin of their healthcare professional
30Gx8mm	Normal weight adults Overweight adults	Overweight adults maybe able to use an 8mm without a lifted skin fold under supervisin of their healthcare professional
29Gx12mm	Overweight adults	Which is a good choice for people with dexterity problems or those who take larger doses. The wider inner diameter makes it easier to press the push button and expel insulin





Safety Insulin Pen Needles

- Visible needle for precision injection.
- Shield sleeve locked immediately after injection.
- Highest level of protection with single-handed technique.
- Wider patient contact surface to reduce pain at injection site.
- Gamma radiation sterilized.
- Fit for most pen injectors.

Red sign shows used and locked



Before used



After used

Needle Length	Color	Packaging
30G x 6mm		100/box ; 30boxes/case
30G x 8mm		100/box ; 30boxes/case
31G x 6mm		100/box ; 30boxes/case
31G x 8mm		100/box ; 30boxes/case



Insulin Syringes

- Specification: 29G 1mL, 0.5mL (12.7mm)
30G 1mL, 0.5mL, 0.3mL (5mm, 8mm)
31G 1mL, 0.5mL, 0.3mL (4mm, 5mm, 8mm)
- Manufactured in class 100k clean rooms.
- ETO sterilized with 5 year shelf life.

Description	Packaging
29G/1mL, 0.5mL (12.7mm)	100/box ; 500/case
30G/1mL, 0.5mL, 0.3mL (5mm)	100/box ; 500/case
30G/1mL, 0.5mL, 0.3mL (8mm)	100/box ; 500/case
31G/1mL, 0.5mL, 0.3mL (4mm)	100/box ; 500/case
31G/1mL, 0.5mL, 0.3mL (5mm)	100/box ; 500/case
31G/1mL, 0.5mL, 0.3mL (8mm)	100/box ; 500/case





Safety Insulin Syringes

- Specification: 29G 1mL, 0.5mL (12.7mm)
30G 1mL, 0.5mL, 0.3mL (5mm, 8mm)
31G 1mL, 0.5mL, 0.3mL (4mm, 5mm, 8mm)
- Sliding sleeve with safety to protect professionals from sharps injuries
- Manufactured in class 100k clean rooms.
- ETO sterilized with 5 year shelf life.

Description	Packaging
29G/1mL, 0.5mL (12.7mm)	100/box ; 500/case
30G/1mL, 0.5mL, 0.3mL (5mm)	100/box ; 500/case
30G/1mL, 0.5mL, 0.3mL (8mm)	100/box ; 500/case
31G/1mL, 0.5mL, 0.3mL (4mm)	100/box ; 500/case
31G/1mL, 0.5mL, 0.3mL (5mm)	100/box ; 500/case
31G/1mL, 0.5mL, 0.3mL (8mm)	100/box ; 500/case

Twist off Lancets

MEDI Lan II Pull Top Lancets

MEDI Lan I



MEDI Lan I+



- Tri bevel sharpened cannula
- Sizes available from 21g to 33g
- Compatible with most lancing devices
- Manufactured in class 100k clean rooms
- Gamma sterilized with 5 year shelf life



- Specification: 28G, 30G, 33G
- Compatible with most lancing devices
- Manufactured in class 100k clean rooms
- Gamma sterilized with 5 year shelf life
- No bent needles when pull off the caps



Carelife USA Inc.
Tel: 770 935 6499
e-mail: info@carelifeusa.com

Shanghai Carelife Int'l Trading Co.,Ltd .
Tel: +86 21 5856 2656
e-mail: sales@shanghaicarelife.com

ReadyLance



Pressure Activated Safety Lancet

ReadyLance Single Use Lancing Device is designed to make blood sampling simple and easy. Touch activation means the lancet can be used in 3 easy steps.

Single use - No need for a lancing device, reduces trans-infection

Safe - Needle retracts immediately after use which reduces the risk of unintentional injuries

Easy - Twist off and press, no need to load lancing device and remove it

Fast - One unit replaces two items which reduces time from loading and unloading the lancets

Accurate - Controlled puncture depth reduces the pain and increases accuracy

ReadyLance Single Use Lancing Device is available in 7 different needle sizes to allow you to choose a suitable lancet or every test you perform.

Color	Gauge/Depth	Amount of bleeding
	30G/1.6mm	
	28G/1.8mm	
	26G/1.8mm	
	23G/2.2mm	
	21G/2.2mm	
	21G/2.2mm Blade	
	18G/2.2mm Blade	



Auto Safety Lancets

- Specification: 21G (2.4mm Blade, 2.0mm Blade, 2.2mm, 1.8mm), 23Gx1.8mm, 26Gx1.8mm, 28Gx1.8mm, 30Gx1.8mm
- Simple one-step function.
- Non-reusable safety feature.
- Custom colors and packing available.
- Gamma sterilized with 5 year shelf life.

Color	Gauge/Depth	Amount of bleeding	Packaging
	30G/1.8mm		100/box ; 1000/case
	28G/1.8mm		100/box ; 1000/case
	26G/1.8mm		100/box ; 1000/case
	23G/1.8mm		100/box ; 1000/case
	21G/1.8mm		100/box ; 1000/case
	21G/2.2mm		100/box ; 1000/case

Mejor




Economic - Safety lancet with simplest structure to provide convenient, reliable, cost-effective blood sampling

Safe - Single use, self-destructing, No possibility of reuse

Easy - intuitive 2 step operation. Blood sampling easier than ever before


Sterile - needles are sterilized with Gamma irradiation

TWIST




Twist off and discard lancet cap. DO NOT PULL.

PRESS



Press platform firmly against sample site. Lancet will automatically retract into body of device.

DISPOSE




Dispose of in a suitable sharps container.

Pressure Activated Safety Lancet

Mejor Single Use Lancing Device is designed to make blood sampling simple and easy. Touch activation means the lancet can be used in 3 easy steps.



- Medium, 2-PLY Prep Pads.
- Isopropyl Alcohol, 70% by Volume.
- For External Antiseptic Use Only.
- Apply topically as needed to cleanse intended area.


Size	Packaging
 2-PLY	1.1 x 2.6 in(2.7 x 6.6cm) 100/box ; 10000/case

Alcohol Prep Pads



Carelife USA Inc.
Tel: 770 935 6499
e-mail: info@carelifeusa.com

Shanghai Carelife Int'l Trading Co.,Ltd .
Tel: +86 21 5856 2656
e-mail: sales@shanghaicarelife.com



Disposable Syringes

- Specification: 1mL, 2mL, 3mL, 5mL, 10mL, 20mL, 30mL, 35mL, 60mL
- Luer lock tip, Luer slip tip or Catheter tip, With or Without needle.
- High clarity barrel, clear scale markings, easy to read.
- Non-toxic and Non-pyrogenic, ETO sterilized with 5 year shelf life.


Description	Packaging
1mL Luer Lock/Slip Syringe without needle	100/box ; 800/case
1mL Luer Lock/Slip Syringe with needle	100/box ; 800/case
2/3mL Luer Lock/Slip Syringe without needle	100/box ; 1200/case
2/3mL Luer Lock/Slip Syringe with needle	100/box ; 1200/case
5/6mL Luer Lock/Slip Syringe without needle	100/box ; 600/case
10mL Luer Lock/Slip Syringe without needle	100/box ; 600/case
10mL Luer Lock/Slip Syringe with needle	100/box ; 600/case
20mL Luer Lock/Slip Syringe without needle	100/box ; 300/case
30mL Luer Lock/Slip Syringe without needle	100/box ; 300/case
35mL Luer Lock/Slip Syringe without needle	100/box ; 300/case
60mL Luer Lock/Slip Syringe without needle	100/box ; 150/case



1mL Luer Lock syringe
Polycarbonate

- 1 mL Luer Lock Syringe is ideal for critical clinical applications where avoiding potential needle disengagement, medication leakage, spray or tubing disconnect is crucial.
- Luer lock hub assures secure connection to needles or IV ports.
- Crystal-clear barrel allows easy measurement of the volume.
- High clarity barrel, clear scale markings, easy to read.

Description	Packaging
1mL PC Syringe	100/box ; 800/case



Safety Needles

Medi®Point Safety Needle is a needle-based safety device that helps you easily convert from non-safety to safety, protecting healthcare workers from exposed to needlestick injuries.



After injection finished, push or slightly knock on the safety cap to protect the needle tip

Description	Packaging
16G x 1" ; 1 1/2"	100/box ; 1000/case
18G x 1" ; 1 1/2"	100/box ; 1000/case
20G x 1" ; 1 1/2"	100/box ; 1000/case
21G x 1" ; 1 1/2" ; 5/8"	100/box ; 1000/case
22G x 1" ; 1 1/2"	100/box ; 1000/case
23G x 1" ; 1 1/2"	100/box ; 1000/case
25G x 1" ; 1 1/2" ; 5/8"	100/box ; 1000/case
27G x 1 1/2" ; 1 1/4"	100/box ; 1000/case
30G x 1/2"	100/box ; 1000/case



Tamper Evident Cap for Syringe

- The pharmacist can feel confident that the integrity of the filled syringe is protected
- Ideal for use with syringes containing controlled substances or chemo drugs
- Unused medication syringes can be safely returned to the pharmacy and re-dispensed, reducing the cost of medication waste
- Supports your pharmacy's efforts to provide a high quality product to the patient

Description	Packaging
Tamper Evident Cap for Syringe	5/bag ; 500/case





Oral Syringes

- Specification: 3mL, 10mL, 12mL, 35mL, 60mL
- Clear barrel allows easy measurement of the volume.
- High clarity barrel, clear scale markings, easy to read.
- Free and smooth movement of plunger.
- Medical graded Polypropylene plunger.



Oral Syringes *Amber*

- Specification: 3mL, 10mL, 12mL, 35mL, 60mL
- Clear barrel allows easy measurement of the volume.
- High clarity barrel, clear scale markings, easy to read.
- Free and smooth movement of plunger.
- Medical graded Polypropylene plunger.



Enteral Syringes

- Specification: 1mL, 3mL, 6mL, 12mL, 20mL, 35mL, 60mL
- Clear barrel allows easy measurement of the volume.
- High clarity barrel, clear scale markings, easy to read.
- Free and smooth movement of plunger.
- Medical graded Polypropylene plunger. Color could be changed.



Enteral Syringes

- Specification: 1mL, 3mL, 6mL, 12mL, 20mL, 35mL, 60mL
- Clear barrel allows easy measurement of the volume.
- High clarity barrel, clear scale markings, easy to read.
- Free and smooth movement of plunger.
- Medical graded Polypropylene plunger. Color could be changed.



Blunt Fill Needles

- Specification: 18G x 1½” , 18G x 2” , 19G x 1½” , 19G x 2”
- Specifically designed for draw-up of medication, reconstitution of medication and accessing IV bags.
- 45° bevel for easy penetration and avoids needles stick injuries.
- Standard Luer hub fits on all size of standard syringes.
- **Blunt Fill Needle With Filter** features a 5 micron membrane designed to prevent large particles from being drawn into the syringe with preparing medications from ampoules.


Description	Packaging
18G x 1 1/2 ”	100/box ; 500/case
18G x 2 ”	100/box ; 500/case
19G x 1 1/2 ”	100/box ; 500/case
19G x 2 ”	100/box ; 500/case



Plastic Blunt Needles

- Plastic blunt needle can be used with all types of split septum IV injection sites and vials for needle-free access, available separately or pre-attached to luer adapters.
- Plastic blunt needle insures rapid withdraw of medication from vial and extremely low medication residual left in vial.
- Dual side ports for thorough flushing of injection sites.
- Clear material for enhanced visualization.

Description	Packaging
Plastic Blunt Needle	100/box ; 1000/case



IV Sets

Hospital IV



6000N 92" 20drop/ml, Vented chamber, Robert clamp, GVS flow regulator, needle free Y-site, rotating luer lock



5500 92", 20drops/ml, vented chamber, Robert clamp, GVS flow regulator, rotating luer lock,



5000 92", 20drops/ml, vented chamber, Robert clamp, GVS flow regulator, injection Y-site, rotating luer lock



2500 18", female luer connect, GVS flow regulator, injection Y-site, rotating luer lock

EMS IV



961001N 96", 10drops/ml, non-vented chamber, roller clamp, needle free Y-site, rotating luer lock

961001 96", 10drops/ml, non-vented chamber, roller clamp, injection Y-site, rotating luer lock



722001 72", 20drops/ml, non vented chamber, roller clamp, injection Y-site, rotating luer lock

726001 72", 60drops/ml, non vented chamber, roller clamp, injection Y-site, rotating luer lock



832003 83" 20drops/ml, non vented chamber, roller clamp, injection Y-site, needle free Y-site, rotating luer lock

836003 83" 60drops/ml, non vented chamber, roller clamp, injection Y-site, needle free Y-site, rotating luer lock



722002 72", 20drops/ml, vented chamber, roller clamp, needle free Y-site, rotating luer lock

726002 72", 60drops/ml, vented chamber, roller clamp, needle free Y-site, rotating luer lock



722003 72", 20drops/ml, vented chamber, slide clamp, roller clamp, needle free Y-site, rotating luer lock

726003 72", 60drops/ml, vented chamber, slide clamp, roller clamp, needle free Y-site, rotating luer lock



Safety Blood Collection Needle Set

- **Safety Blood Collection Needle Set** provides a simple, effective way to collect blood, helps reduce the possibility of accidental needlesticks and features one-handed activation technique.
- Available with a protective cap, reducing the risk of back-end needle stick injuries.
- Smooth and comfortable penetraion, reduces the friction, and tissue damage.

Description	Packaging
Safety Blood Collection Needle Set	50/box ; 1000/case



- Size: 27G, 30G
- Length of needle : 8mm - 38mm
- Sharp tri-bevel point, for maximum comfort.
- Screw-in system: inch type, metric type and universal type.
- Unit package: heat-sealed plastic container or label.

Description	Packaging
27G x 25/32/38mm	100/box ; 500/case
30G x 12/21/25mm	100/box ; 500/case

Dental Needles

D6Uf9L1466Q162



- Specification: 3mL, 6mL, 12mL, 20mL
- Ideal for surgical site Irrigation.
- Compatible with all luer style needles and hubs.
- Disposable , non-sterile.

Description	Packaging
3mL	100/box ; 500/case
6mL	100/box ; 500/case
12mL	100/box ; 500/case
20mL	100/box ; 500/case


Luer-Lock Endo Irrigation Syringes



Available in **4** different types of cannulas



Type I
*Close end,
side-port cannula*



Type II
*Close end,
double-bore cannula*



Type III
*Blunt end,
half-cut cannula*



Type IV
Blunt end cannula

Dental Irrigation Needles

- Size: 27G x 1", 30G x 1"
- Universal luer lock style hub.
- Disposable dental irrigation needle use for canal clean.
- Bendable needle tip designed for irrigation of canals, pockets and fistulas.



Flow Tips
FLOW TIPS

- Size: 18G - 30G
- Length of needle: 13-15mm
- Used for cleaning teeth.
- Pre-bented to the proper angle for placing enchants, sealants cements and other flowable materials.
- All types of flow tips to meet your special.

Description	Packaging
18G(1.2mm)	100/bag ; 50000/case
19G(1.1mm)	100/bag ; 50000/case
20G(0.9mm)	100/bag ; 50000/case
21G(0.8mm)	100/bag ; 50000/case
22G(0.7mm)	100/bag ; 50000/case
23G(0.8mm)	100/bag ; 50000/case
25G(0.5mm)	100/bag ; 50000/case
27G(0.4mm)	100/bag ; 50000/case



- Disposable plastic dental impression syringe for single use.
- Ideal for teeth cleaning, impression material injection, family care after operation and the fluoride treatment.
- OEM orders are welcome!

Description	Packaging
12mL	50/box ; 600/case

Curved Tip Syringes



Shanghai Carelife Int'l Trading Co.,Ltd
#1707 Yingqiao Bldg, No.58 Jinxin Rd,
Pudong, Shanghai 201206, China
Tel: 0086-21-58542656
Email: sales@shanghaicarelife.com

Distributed by:

Carelife USA Inc.
1120 Alpha Dr.
Alpharetta, GA 30004 USA
Tel: 770-935-6499 ; 770-935-6588
Email: info@carelifeusa.com